

兔抗 ATRIP(Phospho-Ser68) 多克隆抗体

中文名称：兔抗 ATRIP(Phospho-Ser68) 多克隆抗体

英文名称： Anti-ATRIP(Phospho-Ser68) rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

抗 原： ATRIP(Phospho-Ser68)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes an essential component of the DNA damage checkpoint. The encoded protein binds to single-stranded DNA coated with replication protein A. The protein also interacts with the ataxia telangiectasia and Rad3 related protein kinase, resulting in its accumulation at intranuclear foci induced by DNA damage.
Applications:	WB
Name of antibody:	ATRIP(Phospho-Ser68)
Immunogen:	Peptide sequence around phosphorylation site of Serine 68 (L-A-S(p)-Q-A) derived from Human ATRIP.
Full name:	ATR interacting protein
Synonyms :	0
SwissProt:	Q8WXE1

WB Predicted band size:	86 kDa
WB Positive control:	NIH/3T3 cells lysates
WB Recommended dilution:	500-1000

