

兔抗 ATRIP 多克隆抗体

中文名称：兔抗 ATRIP 多克隆抗体

英文名称： Anti-ATRIP rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： ATRIP

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an essential component of the DNA damage checkpoint. The encoded protein binds to single-stranded DNA coated with replication protein A. The protein also interacts with the ataxia telangiectasia and Rad3 related protein kinase, resulting in its accumulation at intranuclear foci induced by DNA damage. Multiple transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ATRIP
Immunogen:	Fusion protein of human ATRIP
Full name:	ATR interacting protein
SwissProt:	Q8WXE1
ELISA Recommended dilution	5000-10000

n:	
IHC positive control:	Human thyroid cancer and human tonsil
IHC Recommend dilution:	30-150



