

## DSTN 抗原（重组蛋白）

中文名称：DSTN 抗原（重组蛋白）

英文名称：DSTN Antigen (Recombinant Protein)

别名：ADF; ACTDP; HEL32; bA462D18.2

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	destrin (actin depolymerizing factor)
<b>Synonyms:</b>	ADF; ACTDP; HEL32; bA462D18.2
<b>Swissprot:</b>	P60981
<b>Gene Accession:</b>	BC009477
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The product of this gene belongs to the actin-binding proteins ADF family. This family of proteins is responsible for enhancing the turnover rate of actin in vivo. This gene encodes the actin depolymerizing protein that severs actin filaments (F-actin) and binds to actin monomers (G-actin). Two transcript variants encoding distinct isoforms have been identified for this gene.