

兔抗 ATRN 多克隆抗体

中文名称：兔抗 ATRN 多克隆抗体

英文名称：Anti-ATRN rabbit polyclonal antibody

别名：MGCA; DPPT-L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ATRN

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Multiple transcript variants encoding different isoforms exist for this gene. One of the isoforms is a membrane-bound protein with sequence similarity to the mouse mahogany protein, a receptor involved in controlling obesity. The other two isoforms are secreted proteins involved in the initial immune cell clustering during inflammatory responses that may regulate the chemotactic activity of chemokines.
Applications:	ELISA, IHC
Name of antibody:	ATRN
Immunogen:	Synthetic peptide of human ATRN

Full name:	attractin
Synonyms :	MGCA; DPPT-L
SwissProt:	O75882
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	25-100

