

## 兔抗 CSNK1A1 多克隆抗体

中文名称：兔抗 CSNK1A1 多克隆抗体

英文名称：Anti-CSNK1A1 rabbit polyclonal antibody

别名：CK1; CK1a; CKIa; HLCDGP1; PRO2975; HEL-S-77p

相关类别：一抗

储存：冷冻 (-20℃)

宿主：Rabbit

抗原：CSNK1A1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Casein Kinase I (CK1 or CKI) is the name given to a family of kinases consisting of multiple isoforms ( $\alpha$ , $\alpha'$ , $\beta$ , $\gamma$ 1-3, $\delta$ , and $\epsilon$ ) with a conserved N-terminal kinase domain and a variable C-terminal sequence that determines subcellular localization and regulates enzyme activity.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	CSNK1A1
<b>Immunogen:</b>	Synthetic peptide of human CSNK1A1
<b>Full name:</b>	casein kinase 1, alpha 1
<b>Synonyms :</b>	CK1; CK1a; CKIa; HLCDGP1; PRO2975; HEL-S-77p
<b>SwissProt:</b>	P48729

<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human prostate cancer and Human brain
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	A549 cells
<b>WB Recommended dilution:</b>	200-1000



