

兔抗 CSNK1D 多克隆抗体

中文名称：兔抗 CSNK1D 多克隆抗体

英文名称： Anti-CSNK1D rabbit polyclonal antibody

别名： ASPS; HCKID; FASPS2; CKIdelta

相关类别： 一抗

储存： 冷冻 (-20℃)

抗原： CSNK1D

宿主： Rabbit

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the casein kinase I (CKI) gene family whose members have been implicated in the control of cytoplasmic and nuclear processes, including DNA replication and repair. The encoded protein is highly similar to the mouse and rat CK1 delta homologs. Two transcript variants encoding different isoforms have been found for this gene
Applications:	ELISA, IHC
Name of antibody:	CSNK1D
Immunogen:	Fusion protein of human CSNK1D
Full name:	casein kinase 1, delta

Synonyms :	ASPS; HCKID; FASPS2; CKIdelta
SwissProt:	P48730
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain and human liver cancer
IHC Recommend dilution:	50-200

