

DYDC2 抗原（重组蛋白）

中文名称：DYDC2 抗原（重组蛋白）

英文名称：DYDC2 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

Full name:	DPY30 domain containing 2
Swissprot:	Q96IM9
Gene Accession:	BC007374
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of a family of proteins that contains a DPY30 domain. This gene locus overlaps with a closely related gene on the opposite strand. Alternative splicing results in multiple transcript variants.