

## DYNLL2 抗原（重组蛋白）

中文名称：DYNLL2 抗原（重组蛋白）

英文名称：DYNLL2 Antigen (Recombinant Protein)

别名：Dlc2; DNCL1B; RSPH22

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	dynein, light chain, LC8-type 2
<b>Synonyms:</b>	Dlc2; DNCL1B; RSPH22
<b>Swissprot:</b>	Q96FJ2
<b>Gene Accession:</b>	BC010744
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Dynein light chain 2, cytoplasmic is a protein that in humans is encoded by the DYNLL2 gene. Acts as one of several non-catalytic accessory components of the cytoplasmic dynein 1 complex that are thought to be involved in linking dynein to cargos and to adapter proteins that regulate dynein function. Cytoplasmic dynein 1 acts as a motor for the intracellular retrograde motility of vesicles and organelles along microtubules.