

兔抗 CSNK2A1 多克隆抗体

中文名称：兔抗 CSNK2A1 多克隆抗体

英文名称： Anti-CSNK2A1 rabbit polyclonal antibody

别名： CKII, CK2A1

抗原： CSNK2A1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Casein kinase II is a serine/threonine protein kinase that phosphorylates acidic proteins such as casein. The kinase exists as a tetramer and is composed of an alpha, an alpha-prime, and two beta subunits. The alpha subunits contain the catalytic activity while the beta subunits undergo autophosphorylation. The protein encoded by this gene represents the alpha subunit. While this gene is found on chromosome 20, a related transcribed pseudogene is found on chromosome 11. Three transcript variants encoding two different proteins have been found for this gene.
Applications:	ELISA, WB, IHC

Name of antibody:	CSNK2A1
Immunogen:	Fusion protein of human CSNK2A1
WB Recommended dilution:	200-1000
WB Positive control:	Mouse brain tissue , Hela and Jurkat cells
WB Predicted band size:	45 kDa
IHC Recommend dilution:	25-100
IHC positive control:	Human brain and human cervical cancer
ELISA Recommended dilution:	1000-2000
SwissProt:	P68400
Synonyms :	CKII, CK2A1
Full name:	casein kinase 2, alpha 1 polypeptide



