

## 兔抗 CSRNP3 多克隆抗体

中文名称：兔抗 CSRNP3 多克隆抗体

英文名称：Anti-CSRNP3 rabbit polyclonal antibody

别名：TAIP2; TAIP-2; PPP1R73; FAM130A2

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：CSRNP3

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis By similarity.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CSRNP3
<b>Immunogen:</b>	Synthesized peptide derived from internal of human TAIP-2 .
<b>Full name:</b>	cysteine-serine-rich nuclear protein 3
<b>Synonyms :</b>	TAIP2; TAIP-2; PPP1R73; FAM130A2
<b>SwissProt:</b>	Q8WYN3
<b>WB Predicted band size:</b>	65 kDa

<b>WB Positive control:</b>	COS7 cells lysate
<b>WB Recommended dilution:</b>	500-3000

