

兔抗 CSRP2 多克隆抗体

中文名称：兔抗 CSRP2 多克隆抗体

英文名称： Anti-CSRP2 rabbit polyclonal antibody

别名： cysteine and glycine rich protein 2; CRP2; LMO5; SmLIM

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： CSRP2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	CSRP2 is a member of the CSRP family of genes, encoding a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. CRP2 contains two copies of the cysteine-rich amino acid sequence motif (LIM) with putative zinc-binding activity, and may be involved in regulating ordered cell growth. Other genes in the family include CSRP1 and CSRP3. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, WB, IHC
Name of antibody:	CSRP2

Immunogen:	Fusion protein of human CSRP2
Full name:	cysteine and glycine rich protein 2
Synonyms:	CRP2; LMO5; SmLIM
SwissProt:	Q16527
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil and Human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	21 kDa
WB Positive control:	HepG2 and A549 cell lysates
WB Recommended dilution:	1000-5000



