

兔抗 AVPR1B 多克隆抗体

中文名称：兔抗 AVPR1B 多克隆抗体

英文名称： Anti-AVPR1B rabbit polyclonal antibody

别名： arginine vasopressin receptor 1B; AVPR3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： AVPR1B

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene acts as receptor for arginine vasopressin. This receptor belongs to the subfamily of G-protein coupled receptors which includes AVPR1A, V2R and OXT receptors. Its activity is mediated by G proteins which stimulate a phosphatidylinositol-calcium second messenger system. The receptor is primarily located in the anterior pituitary, where it stimulates ACTH release. It is expressed at high levels in ACTH-secreting pituitary adenomas as well as in bronchial carcinoids responsible for the ectopic ACTH syndrome. A spliced antisense transcript of this gene has been reported but its function is not known.

	own.
Applications:	ELISA, IHC
Name of antibody:	AVPR1B
Immunogen:	Synthetic peptide of human AVPR1B
Full name:	arginine vasopressin receptor 1B
Synonyms:	AVPR3
SwissProt:	P47901
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	10-50

