

## 兔抗 AVPR2 多克隆抗体

中文名称：兔抗 AVPR2 多克隆抗体

英文名称：Anti-AVPR2 rabbit polyclonal antibody

别名：DI1, DIR, NDI, V2R, ADHR, DIR3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：AVPR2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the ASPP (apoptosis-stimulating protein of p53) family of p53 interacting proteins. The protein contains four ankyrin repeats and an SH3 domain involved in protein-protein interactions. It is localized to the perinuclear region of the cytoplasm, and regulates apoptosis and cell growth through interactions with other regulatory molecules including members of the p53 family. Multiple transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	AVPR2

<b>Immunogen:</b>	Synthetic peptide of human AVPR2
<b>Full name:</b>	arginine vasopressin receptor 2
<b>Synonyms :</b>	DI1, DIR, NDI, V2R, ADHR, DIR3
<b>SwissProt:</b>	P30518
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human colon cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	25-100

