

EMC8 抗原（重组蛋白）

中文名称：EMC8 抗原（重组蛋白）

英文名称：EMC8 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

| | |
|---------------------------|---|
| Full name: | ER membrane protein complex subunit 8 |
| Synonyms: | NOC4; COX4NB; C16orf2; C16orf4; FAM158B |
| Swissprot: | O43402 |
| Gene Accession: | BC001472 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | COX4NB (Neighbor of COX4) is a 210 amino acid protein encoded by the human gene COX4NB. COX4NB belongs to the UPF0172 (NOC4) family and is found on chromosome 16, adjacent to the gene that encodes COX4. The GAN gene is located on chromosome 16 and, with mutation, may lead to giant axonal neuropathy, a nervous system disorder characterized by increasing malfunction with growth. The rare disorder Rubinstein-Taybi syndrome is also associated with chromosome 16 through the CREBBP gene which encodes a critical CREB binding protein. Crohn's disease is a gastrointestinal inflammatory condition associated with chromosome 16 through |

gh the NOD2 gene. An association with systemic lupus erythemato
sis and a number of other auto-immune disorders with the perice
ntromeric region of chromosome 16 has led to the identification o
f SLC5A11 as a potential autoimmune modifier.