

兔抗 B3GALT6 多克隆抗体

中文名称：兔抗 B3GALT6 多克隆抗体

英文名称： Anti-B3GALT6 rabbit polyclonal antibody

别名： beta-1,3-galactosyltransferase 6; EDSP2; SEMDJL1; beta3GalT6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： B3GALT6

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The enzyme encoded by this intronless gene is a beta-1,3-galactosyltransferase found in the medial Golgi apparatus , where it catalyzes the transfer of galactose from UDP-galactose to substrates containing a terminal beta-linked galactose moiety. The encoded enzyme has a particular affinity for galactose-beta-1,4-xylose found in the linker region of glycosamines. This enzyme is required for glycosaminoglycan synthesis.
Applications:	ELISA, IHC
Name of antibody:	B3GALT6
Immunogen:	Synthetic peptide of human B3GALT6

Full name:	beta-1,3-galactosyltransferase 6
Synonyms:	EDSP2; SEMDJL1; beta3GalT6
SwissProt:	Q96L58
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

