

FHL2 抗原（重组蛋白）

中文名称：FHL2 抗原（重组蛋白）

英文名称：FHL2 Antigen (Recombinant Protein)

别名：DRAL; AAG11; FHL-2; SLIM3; SLIM-3

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	four and a half LIM domains 2
Synonyms:	DRAL; AAG11; FHL-2; SLIM3; SLIM-3
Swissprot:	Q14192
Gene Accession:	NP_963849
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the four-and-a-half-LIM-only protein family. Family members contain two highly conserved, tandemly arranged, zinc finger domains with four highly conserved cysteines binding a zinc atom in each zinc finger. This protein is thought to have a role in the assembly of extracellular membranes. Also, this gene is down-regulated during transformation of normal myoblasts to rhabdomyosarcoma cells and the encoded protein may function as a

link between presenilin-2 and an intracellular signaling pathway. Multiple alternatively spliced variants, encoding the same protein, have been identified.