

FSD1 抗原（重组蛋白）

中文名称：FSD1 抗原（重组蛋白）

英文名称： FSD1 Antigen (Recombinant Protein)

别名： fibronectin type III and SPRY domain containing 1; MIR1; GLFND

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 297-496 amino acids of human FSD1

技术规格

Full name:	fibronectin type III and SPRY domain containing 1
Synonyms:	MIR1; GLFND
Swissprot:	Q9BTV5
Gene Accession:	BC016442
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a centrosome associated protein that is characterized by an N-terminal coiled-coil region downstream of B-box (BBC) domain, a central fibronectin type III domain, and a C-terminal repeats in splA and RyR (SPRY) domain. The encoded protein associates with a subset of microtubules and may be involved in the stability and organization of microtubules during cytokinesis.