

兔抗 BAALC 多克隆抗体

中文名称: 兔抗 BAALC 多克隆抗体

英文名称: Anti-BAALC rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: BAALC

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene was identified by gene expression studies in patients with acute myeloid leukemia (AML). The gene is conserved among mammals and is not found in lower organisms. Tissues that express this gene develop from the neuroectoderm. Multiple alternatively spliced transcript variants that encode different proteins have been described for this gene; however, some of the transcript variants are found only in AML cell lines.
Applications:	ELISA, IHC
Name of antibody:	BAALC
Immunogen:	Full length fusion protein
Full name:	brain and acute leukemia, cytoplasmic
SwissProt:	Q8WXS3

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

