

兔抗 CST9L 多克隆抗体

中文名称：兔抗 CST9L 多克隆抗体

英文名称： Anti-CST9L rabbit polyclonal antibody

别名： CTES7B; bA218C14.1

抗原： CST9L

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. There are three inhibitory families in the superfamily, including the type 1 cystatins (stefins), type 2 cystatins and the kininogens.
Applications:	IHC
Name of antibody:	CST9L
Immunogen:	Synthesized peptide derived from internal of human CST9L.
Full name:	cystatin 9-like
Synonyms :	CTES7B; bA218C14.1

SwissProt:	Q9H4G1
IHC positive control:	Human placenta tissue
IHC Recommend dilution:	50-100

