

## 兔抗 CST11 多克隆抗体

中文名称：兔抗 CST11 多克隆抗体

英文名称： Anti-CST11 rabbit polyclonal antibody

别名： SC13; CST8L; CTES2; dJ322G13.6

抗原： CST11

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. There are three inhibitory families in the superfamily, including the type 1 cystatins (stefins), type 2 cystatins and the kininogens. The type 2 cystatin proteins are a class of cysteine proteinase inhibitors found in a variety of human fluids and secretions.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CST11

<b>Immunogen:</b>	Synthesized peptide derived from internal of human CST1 1.
<b>Full name:</b>	cystatin 11
<b>Synonyms :</b>	SC13; CST8L; CTES2; dJ322G13.6
<b>SwissProt:</b>	Q9H112
<b>WB Predicted band size:</b>	17 kDa
<b>WB Positive control:</b>	A549 cells lysate
<b>WB Recommended dilution:</b>	500-3000

