

## 兔抗 Bad (Phospho-Ser112) 多克隆抗体

中文名称：兔抗 Bad (Phospho-Ser112) 多克隆抗体

英文名称：Anti-Bad (Phospho-Ser112) rabbit polyclonal antibody

别名：Bbc2; A1325008

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：Bad (Phospho-Ser112)

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

The protein encoded by BAD gene is a member of the BCL-2 family. BCL-2 family members are known to be regulators of programmed cell death. This protein positively regulates cell apoptosis by forming heterodimers with BCL-xL and BCL-2, and reversing their death repressor activity. Proapoptotic activity of this protein is regulated through its phosphorylation. Protein kinases AKT and MAP kinase, as well as protein phosphatase calcineurin were found to be involved in the regulation of this protein. Alternative splicing of this gene results in two transcript variants which encode the same isoform.

<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	Bad (Phospho-Ser112)
<b>Immunogen:</b>	Synthetic peptide of mouse Bad (Phospho-Ser112)
<b>Full name:</b>	BCL2-associated agonist of cell death (Phospho-Ser112)
<b>Synonyms :</b>	Bbc2; AI325008
<b>SwissProt:</b>	Q61337
<b>IHC positive control:</b>	Human breast carcinoma
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	Cos7 cells treated with EGF
<b>WB Recommended dilution:</b>	500-1000



