

## 兔抗 CSTF1 多克隆抗体

中文名称：兔抗 CSTF1 多克隆抗体

英文名称： Anti-CSTF1 rabbit polyclonal antibody

别名： CstF-50; CstFp50

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： CSTF1

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes one of three subunits which combine to form cleavage stimulation factor (CSTF). CSTF is involved in the polyadenylation and 3'end cleavage of pre-mRNAs. Similar to mammalian G protein beta subunits, this protein contains transducin-like repeats. Several transcript variants with different 5' UTR, but encoding the same protein, have been found for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	CSTF1
<b>Immunogen:</b>	Synthesized peptide derived from N-terminal of huma

	n CSTF1.
<b>Full name:</b>	cleavage stimulation factor, 3' pre-RNA, subunit 1, 50 kDa
<b>Synonyms :</b>	CstF-50; CstFp50
<b>SwissProt:</b>	Q05048
<b>WB Predicted band size:</b>	48 kDa
<b>WB Positive control:</b>	293 cells and Huvec cells lysates
<b>WB Recommended dilution:</b>	500-3000

