

兔抗 Bad(Phospho-Ser155) 多克隆抗体

- 中文名称： 兔抗 Bad(Phospho-Ser155) 多克隆抗体
- 英文名称： Anti-Bad(Phospho-Ser155) rabbit polyclonal antibody
- 别名： Bbc2; A1325008
- 相关类别： 一抗
- 储存： 冷冻（-20℃） 避光
- 宿主： Rabbit
- 抗原： Bad(Phospho-Ser155)
- 反应种属： Human Mouse Rat
- 标记物： Unconjugate
- 克隆类型： Unconjugate

技术规格

Background:	Promotes cell death. Successfully competes for the binding to Bcl-X(L), Bcl-2 and Bcl-W, thereby affecting the level of heterodimerization of these proteins with BAX. Can reverse the death repressor activity of Bcl-X(L), but not that of Bcl-2. Appears to act as a link between growth factor receptor signaling and the apoptotic pathways.
Applications:	WB, IHC
Name of antibody:	Bad(Phospho-Ser155)
Immunogen:	Peptide sequence around phosphorylation site of serine 155(R-M-S(p)-D-E) derived from Human BAD.
Full name:	BCL2-associated agonist of cell death

Synonyms :	Bbc2; AI325008
SwissProt:	Q61337
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	22 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution:	500-1000



