

## GNG2 抗原（重组蛋白）

中文名称：GNG2 抗原（重组蛋白）

英文名称：GNG2 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	guanine nucleotide binding protein (G protein), gamma 2
<b>Swissprot:</b>	P59768
<b>Gene Accession:</b>	BC020774
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes one of the gamma subunits of a guanine nucleotide-binding protein. Such proteins are involved in signaling mechanisms across membranes. Various subunits forms heterodimers which then interact with the different signal molecules.