

GPN2 抗原（重组蛋白）

中文名称：GPN2 抗原（重组蛋白）

英文名称：GPN2 Antigen (Recombinant Protein)

别名：GPN-loop GTPase 2; ATPBD1B

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 2-310 amino acids of human GPN2

技术规格

Full name:	GPN-loop GTPase 2
Synonyms:	ATPBD1B
Swissprot:	Q9H9Y4
Gene Accession:	BC008634
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Small GTPase required for proper localization of RNA polymerase II and III (RNAPII and RNAPIII). May act at an RNAP assembly step prior to nuclear import.