

兔抗 CSTL1 多克隆抗体

中文名称：兔抗 CSTL1 多克隆抗体

英文名称： Anti-CSTL1 rabbit polyclonal antibody

别名： CTES1; RCET11; dJ322G13.4

抗原： CSTL1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The cystatin superfamily encompasses proteins that contain multiple cystatin-like sequences. Some of the members are active cysteine protease inhibitors, while others have lost or perhaps never acquired this inhibitory activity. There are three inhibitory families in the superfamily, including the type 1 cystatins (stefins), type 2 cystatins and the kininogens. The type 2 cystatin proteins are a class of cysteine proteinase inhibitors found in a variety of human fluids and secretions.
Applications:	IF
Name of antibody:	CSTL1
Immunogen:	Synthesized peptide derived from internal of human CSTL1.

Full name:	cystatin-like 1
Synonyms :	CTES1; RCET11; dJ322G13.4
SwissProt:	Q9H114
IF positive control:	A549 cells
IF Recommend dilution:	100-500

