

GRP 抗原(重组蛋白)

- 中文名称: GRP 抗原(重组蛋白)
- 英文名称: GRP Antigen (Recombinant Protein)
- 储存: 冷冻(-20℃)
- 相关类别: 抗原

概述

Full length fusion protein

技术规格

Full name:	gastrin-releasing peptide
Synonyms:	BN; GRP-10; proGRP; preproGRP
Swissprot:	P07492
Gene Accession:	BC004488
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the bombesin-like family of gastrin -releasing peptides. Its preproprotein, following cleavage of a signal peptide, is further processed to produce either the 27 aa gastrin-rel easing peptide or the 10 aa neuromedin C. These smaller peptides r egulate numerous functions of the gastrointestinal and central nervo us systems, including release of gastrointestinal hormones, smooth muscle cell contraction, and epithelial cell proliferation. These peptid es are also likely to play a role in human cancers of the lung, colon , stomach, pancreas, breast, and prostate. Alternative splicing results in multiple transcript variants encoding different isoforms.

