

兔抗 CTBP1(Ab-422) 多克隆抗体

中文名称：兔抗 CTBP1(Ab-422) 多克隆抗体

英文名称：Anti-CTBP1(Ab-422) rabbit polyclonal antibody

别名：BARS

相关类别：一抗

抗原：CTBP1(Ab-422)

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Corepressor targeting diverse transcription regulators such as GLIS2 or BCL6. Has dehydrogenase activity. Involved in controlling the equilibrium between tubular and stacked structures in the Golgi complex. Functions in brown adipose tissue (BAT) differentiation.
Applications:	WB, IHC
Name of antibody:	CTBP1(Ab-422)
Immunogen:	Synthesized non-phosphopeptide derived from human CtBP1 around the phosphorylation site of serine 422 (A-P-S(p)-P-G).
Full name:	C-terminal binding protein 1

Synonyms :	BARS
SwissProt:	Q13363
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	48 kDa
WB Positive control:	Jurkatcells lysate
WB Recommended dilution:	500-3000



