

## 兔抗 BAG6 多克隆抗体

中文名称：兔抗 BAG6 多克隆抗体

英文名称：Anti-BAG6 rabbit polyclonal antibody

别名：G3; BAT3; BAG-6; D6S52E

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BAG6

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Background:</b>       | This gene was first characterized as part of a cluster of genes located within the human major histocompatibility complex class III region. This gene encodes a nuclear protein that is cleaved by caspase 3 and is implicated in the control of apoptosis. In addition, the protein forms a complex with E1A binding protein p300 and is required for the acetylation of p53 in response to DNA damage. Multiple transcript variants encoding different isoforms have been found for this gene. |
| <b>Applications:</b>     | ELISA, WB                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Name of antibody:</b> | BAG6                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |

|                                    |                                 |
|------------------------------------|---------------------------------|
| <b>Immunogen:</b>                  | Synthetic peptide of human BAG6 |
| <b>Full name:</b>                  | BCL2-associated athanogene 6    |
| <b>Synonyms :</b>                  | G3; BAT3; BAG-6; D6S52E         |
| <b>SwissProt:</b>                  | P46379                          |
| <b>ELISA Recommended dilution:</b> | 1000-2000                       |
| <b>WB Predicted band size:</b>     | 119 kDa                         |
| <b>WB Positive control:</b>        | 293T, hela and hepG2 cells      |
| <b>WB Recommended dilution:</b>    | 200-1000                        |

