

兔抗 CTBP2 多克隆抗体

中文名称：兔抗 CTBP2 多克隆抗体

英文名称： Anti-CTBP2 rabbit polyclonal antibody

抗 原： CTBP2

储 存： 冷冻（-20℃）

宿 主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene produces alternative transcripts encoding two distinct proteins. One protein is a transcriptional repressor, while the other isoform is a major component of specialized synapses known as synaptic ribbons. Both proteins contain a NAD ⁺ binding domain similar to NAD ⁺ -dependent 2-hydroxyacid dehydrogenases. A portion of the 3' untranslated region was used to map this gene to chromosome 21q21.3; however, it was noted that similar loci elsewhere in the genome are likely. Blast analysis shows that this gene is present on chromosome 10.
Applications:	ELISA, IHC
Name of antibody:	CTBP2
Immunogen:	Fusion protein of human CTBP2

Full name:	C-terminal binding protein 2
SwissProt:	P56545
IHC positive control:	Human gastric cancer and human ovarian cancer
IHC Recommend dilution:	25-100

