

兔抗 BAIAP2L2 多克隆抗体

中文名称：兔抗 BAIAP2L2 多克隆抗体

英文名称： Anti-BAIAP2L2 rabbit polyclonal antibody

抗 原： BAIAP2L2

储 存： 冷冻（-20℃） 避光

宿 主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene binds phosphoinositides and promotes the formation of planar or curved membrane structures. The encoded protein is found in RAB13-positive vesicles and at intercellular contacts with the plasma membrane.
Applications:	WB
Name of antibody:	BAIAP2L2
Immunogen:	Synthesized peptide derived from internal of human BAIAP2L2.
Full name:	BAI1-associated protein 2-like 2
Synonyms :	0
SwissProt:	Q6UXY1
WB Predicted band size:	59 kDa

WB Positive control:	RAW264.7 cells lysate
WB Recommended dilution:	500-3000

