

兔抗 BAIAP2 多克隆抗体

中文名称：兔抗 BAIAP2 多克隆抗体

英文名称： Anti-BAIAP2 rabbit polyclonal antibody

别名： BAP2, FLAF3, IRSP53

抗原： BAIAP2

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The protein encoded by this gene has been identified as a brain-specific angiogenesis inhibitor (BAI1)-binding protein. This adaptor protein links membrane bound G-proteins to cytoplasmic effector proteins. This protein functions as an insulin receptor tyrosine kinase substrate and suggests a role for insulin in the central nervous system. It also associates with a downstream effector of Rho small G proteins, which is associated with the formation of stress fibers and cytokinesis. This protein is involved in lamellipodia and filopodia formation in motile cells and may affect neuronal growth-cone guidance. This protei

	n has also been identified as interacting with the dentat orubral-pallidolusian atrophy gene, which is associated with an autosomal dominant neurodegenerative disease. Alternative splicing results in multiple transcript variants encoding distinct isoforms
Applications:	ELISA, IHC
Name of antibody:	BAIAP2
Immunogen:	Fusion protein of human BAIAP2
Full name:	BAI1-associated protein 2
Synonyms :	BAP2, FLAF3, IRSP53
SwissProt:	Q9UQB8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer
IHC Recommend dilution:	50-200

