

## 兔抗 BAK1 多克隆抗体

中文名称：兔抗 BAK1 多克隆抗体

英文名称： Anti-BAK1 rabbit polyclonal antibody

别名： BAK; CDN1; BCL2L7; BAK-LIKE

抗原： BAK1

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型 rabbit polyclonal

### 技术规格

**Background:**

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form oligomers or heterodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. This protein localizes to mitochondria, and functions to induce apoptosis. It interacts with and accelerates the opening of the mitochondrial voltage-dependent anion channel, which leads to a loss in membrane potential and the release of cytochrome c. This protein also interacts with the tumor suppressor P53 after exposure to cell stress.

<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	BAK1
<b>Immunogen:</b>	Fusion protein of human BAK1
<b>Full name:</b>	BCL2-antagonist/killer 1
<b>Synonyms :</b>	BAK; CDN1; BCL2L7; BAK-LIKE
<b>SwissProt:</b>	Q16611
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human esophagus cancer and human breast cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	23 kDa
<b>WB Positive control:</b>	Hela, A431, 293T and Jurkat cells
<b>WB Recommended dilution:</b>	1000-5000



