

## 兔抗 BARX1 多克隆抗体

中文名称：兔抗 BARX1 多克隆抗体

英文名称：Anti-BARX1 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：BARX1

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the Bar subclass of homeobox transcription factors. Studies of the mouse and chick homolog suggest the encoded protein may play a role in developing teeth and craniofacial mesenchyme of neural crest origin. The protein may also be associated with differentiation of stomach epithelia.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	BARX1
<b>Immunogen:</b>	Fusion protein of human BARX1
<b>Full name:</b>	BARX homeobox 1
<b>SwissProt:</b>	Q9HBU1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer

<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	27 kDa
<b>WB Positive control:</b>	Mouse stomach tissue lysate
<b>WB Recommended dilution:</b>	200-1000



