

兔抗 CTNNB1(Phospho-Tyr489) 多克隆抗体

中文名称：兔抗 CTNNB1(Phospho-Tyr489) 多克隆抗体

英文名称： Anti-CTNNB1(Phospho-Tyr489) rabbit polyclonal antibody

别名： CTNNB; MRD19; armadillo

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： CTNNB1(Phospho-Tyr489)

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The protein encoded by this gene is part of a complex of proteins that constitute adherens junctions (AJs). AJs are necessary for the creation and maintenance of epithelial cell layers by regulating cell growth and adhesion between cells. The encoded protein also anchors the actin cytoskeleton and may be responsible for transmitting the contact inhibition signal that causes cells to stop dividing once the epithelial sheet is complete.
Applications:	WB
Name of antibody:	CTNNB1(Phospho-Tyr489)
Immunogen:	Peptide sequence around phosphorylation site of Tyrosi

	ne 489 (L-H-Y(p)-G-L) derived from Human Catenin- beta a.
Full name:	catenin (cadherin-associated protein), beta 1, 88kDa
Synonyms :	CTNNB; MRD19; armadillo
SwissProt:	P35222
WB Predicted band size:	85 kDa
WB Positive control:	Mouse renin lysates
WB Recommended dilution:	500-1000

