

## 兔抗 CTNNBIP1 多克隆抗体

中文名称：兔抗 CTNNBIP1 多克隆抗体

英文名称：Anti-CTNNBIP1 rabbit polyclonal antibody

别名：ICAT

相关类别：一抗

储存：冷冻（-20℃）

抗原：CTNNBIP1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene binds CTNNB1 and prevents interaction between CTNNB1 and TCF family members. The encoded protein is a negative regulator of the Wnt signaling pathway. Two transcript variants encoding the same protein have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CTNNBIP1
<b>Immunogen:</b>	Fusion protein of human CTNNBIP1
<b>Full name:</b>	catenin, beta interacting protein 1
<b>Synonyms :</b>	ICAT
<b>SwissProt:</b>	Q9NSA3

<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human breast cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200

