

HMGN3 抗原（重组蛋白）

中文名称：HMGN3 抗原（重组蛋白）

英文名称：HMGN3 Antigen (Recombinant Protein)

别名：TRIP7; PNAS-24; PNAS-25

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	high mobility group nucleosomal binding domain 3
Synonyms:	TRIP7; PNAS-24; PNAS-25
Swissprot:	Q15651
Gene Accession:	BC009529
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	Thyroid hormone receptors are hormone-dependent transcription factors that regulate expression of a variety of specific target genes. The protein encoded by this gene binds thyroid hormone receptor beta, but only in the presence of thyroid hormone. The encoded protein, a member of the HMGN protein family, is thought to reduce the compactness of the chromatin fiber in nucleosomes, thereby enhancing transcription from chromatin templates.