

## 兔抗 CTNND1(Ab-228) 多克隆抗体

中文名称：兔抗 CTNND1(Ab-228) 多克隆抗体

英文名称： Anti-CTNND1(Ab-228) rabbit polyclonal antibody

别名： CAS; p120; CTNND; P120CAS; P120CTN; p120(CAS); p120(CTN)

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： CTNND1(Ab-228)

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the Armadillo protein family, which function in adhesion between cells and signal transduction. Multiple translation initiation codons and alternative splicing result in many different isoforms being translated. Not all of the full-length nature of the described transcript variants have been determined. Read-through transcription also exists between this gene and the neighboring upstream thioredoxin-related transmembrane protein 2 (TMX2) gene.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CTNND1(Ab-228)

<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human Catenin- $\delta$ 1 around the phosphorylation site of tyrosine 228 (D-N-Y(p)-G-S).
<b>Full name:</b>	catenin (cadherin-associated protein), delta 1
<b>Synonyms :</b>	CAS; p120; CTNND; P120CAS; P120CTN; p120(CAS); p120(CTN)
<b>SwissProt:</b>	O60716
<b>IHC positive control:</b>	Human brain tissue
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	108 kDa
<b>WB Positive control:</b>	Huvec cells lysate
<b>WB Recommended dilution:</b>	500-3000



