

HOXA6 抗原（重组蛋白）

中文名称：HOXA6 抗原（重组蛋白）

英文名称：HOXA6 Antigen (Recombinant Protein)

相关类别：抗原

储存：冷冻（-20℃）

概述

Full length fusion protein

技术规格

Full name:	homeobox A6
Synonyms:	HOX1, HOX1B, HOX1.2
Swissprot:	P31267
Gene Accession:	BC069497
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	In vertebrates, the genes encoding the class of transcription factors called homeobox genes are found in clusters named A, B, C, and D on four separate chromosomes. Expression of these proteins is spatially and temporally regulated during embryonic development. This gene is part of the A cluster on chromosome 7 and encodes a DNA-binding transcription factor which may regulate gene expression, morphogenesis, and differentiation.