

## HOXB1 抗原（重组蛋白）

中文名称：HOXB1 抗原（重组蛋白）

英文名称：HOXB1 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

<b>Full name:</b>	homeobox B1
<b>Synonyms:</b>	HOX2; HCFP3; HOX2I; Hox-2.9
<b>Swissprot:</b>	P14653
<b>Gene Accession:</b>	BC096192
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, located on different chromosomes, consisting of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXB genes located in a cluster on chromosome 17.