

HSF2 抗原（重组蛋白）

中文名称：HSF2 抗原（重组蛋白）

英文名称：HSF2 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

| | |
|---------------------------|--|
| Full name: | heat shock transcription factor 2 |
| Synonyms: | HSF 2, HSTF 2 |
| Swissprot: | Q03933 |
| Gene Accession: | BC005329 |
| Purity: | >85%, as determined by Coomassie blue stained SDS-PAGE |
| Expression system: | Escherichia coli |
| Tags: | His tag C-Terminus, GST tag N-Terminus |
| Background: | The protein encoded by this gene belongs to the HSF family of transcription factors that bind specifically to the heat-shock promoter element and activate transcription. Heat shock transcription factors activate heat-shock response genes under conditions of heat or other stresses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. |