

兔抗 CTSS 多克隆抗体

中文名称: 兔抗 CTSS 多克隆抗体

英文名称: Anti-CTSS rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: CTSS

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene, a member of the peptidase C1 family, is a lysosomal cysteine proteinase that may participate in the degradation of antigenic proteins to peptides for presentation on MHC class II molecules. The encoded protein can function as an elastase over a broad pH range in alveolar macrophages. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	CTSS
Immunogen:	Synthetic peptide of human CTSS
Full name:	cathepsin S
SwissProt:	P25774

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-100
WB Predicted band size:	37 kDa
WB Positive control:	Jurkat, 231 cell, Mouse brain tissue lysates
WB Recommended dilution:	500-2000



