

## 兔抗 BCAN 多克隆抗体

中文名称：兔抗 BCAN 多克隆抗体

英文名称：Anti-BCAN rabbit polyclonal antibody

别名：brevican; BEHAB; CSPG7

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BCAN

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the lectican family of chondroitin sulfate proteoglycans that is specifically expressed in the central nervous system. This protein is developmentally regulated and may function in the formation of the brain extracellular matrix. This protein is highly expressed in gliomas and may promote the growth and cell motility of brain tumor cells. Alternate splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BCAN
<b>Immunogen:</b>	Fusion protein of human BCAN

<b>Full name:</b>	brevican
<b>Synonyms:</b>	BEHAB; CSPG7
<b>SwissProt:</b>	Q96GW7
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and human gastric cancer
<b>IHC Recommend dilution:</b>	30-150



