

IL3 抗原（重组蛋白）

中文名称：IL3 抗原（重组蛋白）

英文名称：IL3 Antigen (Recombinant Protein)

别名：IL-3; MCGF; MULTI-CSF

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	Interleukin 3 (colony-stimulating factor, multiple)
Synonyms:	IL-3; MCGF; MULTI-CSF
Swissprot:	P08700
Gene Accession:	BC069472
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The protein encoded by this gene is a potent growth promoting cytokine. This cytokine is capable of supporting the proliferation of a broad range of hematopoietic cell types. It is involved in a variety of cell activities such as cell growth, differentiation and apoptosis. This cytokine has been shown to also possess neurotrophic activity, and it may be associated with neurologic disorders.