

ING2 抗原（重组蛋白）

中文名称：ING2 抗原（重组蛋白）

英文名称：ING2 Antigen (Recombinant Protein)

别名：ING1L; p33ING2

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	inhibitor of growth family, member 2
Synonyms:	ING1L; p33ING2
Swissprot:	Q9H160
Gene Accession:	BC030128
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene is a member of the inhibitor of growth (ING) family. Members of the ING family associate with and modulate the activity of histone acetyltransferase (HAT) and histone deacetylase (HDAC) complexes and function in DNA repair and apoptosis.