

IRGM 抗原（重组蛋白）

中文名称：IRGM 抗原（重组蛋白）

英文名称：IRGM Antigen (Recombinant Protein)

别名：IFI1; IRGM1; LRG47; LRG-47

储存：冷冻（-20℃）

相关类别：抗原

概述

Full length fusion protein

技术规格

Full name:	immunity-related GTPase family, M
Synonyms:	IFI1; IRGM1; LRG47; LRG-47
Swissprot:	A1A4Y4
Gene Accession:	BC128168
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the p47 immunity-related GTPase family. The encoded protein may play a role in the innate immune response by regulating autophagy formation in response to intracellular pathogens. Polymorphisms that affect the normal expression of this gene are associated with a susceptibility to Crohn's disease and tuberculosis.