

兔抗 BCAP31 多克隆抗体

中文名称：兔抗 BCAP31 多克隆抗体

英文名称： Anti-BCAP31 rabbit polyclonal antibody

抗 原： BCAP31

别 名： CDM; BAP31; 6C6-AG; DXS1357E

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene encodes a member of the B-cell receptor associated protein 31 superfamily. The encoded protein is a multi-pass transmembrane protein of the endoplasmic reticulum that is involved in the anterograde transport of membrane proteins from the endoplasmic reticulum to the Golgi and in caspase 8-mediated apoptosis. Microdeletions in this gene are associated with contiguous ABCD1/DXS1375E deletion syndrome (CADD5), a neonatal disorder. Alternative splicing of this gene results in multiple transcript variants. Two related pseudogenes have been identified on chromosome 16.

Applications:	ELISA, WB, IHC
Name of antibody:	BCAP31
Immunogen:	Fusion protein of human BCAP31
Full name:	B-cell receptor-associated protein 31
Synonyms:	CDM; BAP31; 6C6-AG; DXS1357E
SwissProt:	P51572
IHC positive control:	Human gastric cancer and human thyroid cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	28 kDa
WB Positive control:	Human placenta tissue and A549 cells, Raji and hela c ells
WB Recommended dilution:	500-2000



