

兔抗 CTSZ 多克隆抗体

中文名称：兔抗 CTSZ 多克隆抗体

英文名称： Anti-CTSZ rabbit polyclonal antibody

别名： CTSX

抗原： CTSZ

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a lysosomal cysteine proteinase and member of the peptidase C1 family. It exhibits both carboxy-monopeptidase and carboxy-dipeptidase activities. The encoded protein has also been known as cathepsin X and cathepsin P. This gene is expressed ubiquitously in cancer cell lines and primary tumors and, like other members of this family, may be involved in tumorigenesis.
Applications:	ELISA, WB, IHC
Name of antibody:	CTSZ
Immunogen:	Fusion protein of human CTSZ

Full name:	Cathepsin Z
Synonyms :	CTSX
SwissProt:	Q9UBR2
ELISA Recommended dilution:	500-2000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	10-50
WB Predicted band size:	34 kDa
WB Positive control:	Hela cells
WB Recommended dilution:	500-1000



