

兔抗 BCHE 多克隆抗体

中文名称：兔抗 BCHE 多克隆抗体

英文名称：Anti-BCHE rabbit polyclonal antibody

别名：butyrylcholinesterase; E1; CHE1; CHE2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BCHE

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Mutant alleles at the BCHE locus are responsible for suxamethonium sensitivity. Homozygous persons sustain prolonged apnea after administration of the muscle relaxant suxamethonium in connection with surgical anesthesia. The activity of pseudocholinesterase in the serum is low and its substrate behavior is atypical. In the absence of the relaxant, the homozygote is at no known disadvantage.
Applications:	ELISA, IHC
Name of antibody:	BCHE
Immunogen:	Fusion protein of human BCHE

Full name:	butyrylcholinesterase
Synonyms:	E1; CHE1; CHE2
SwissProt:	P06276
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human thyroid cancer and human gastric cancer
IHC Recommend dilution:	40-200



