

## 兔抗 CUEDC2 多克隆抗体

中文名称：兔抗 CUEDC2 多克隆抗体

英文名称：Anti-CUEDC2 rabbit polyclonal antibody

别名：C10orf66; bA18I14.5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：CUEDC2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	CUEDC2 is involved in ubiquitin- and proteasome-mediated degradation of progesterone receptor and estrogen receptor (ER)-alpha. CUEDC2 reduced PRB protein content and promoted progesterone-induced PRB degradation via the ubiquitin-proteasome pathway.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	CUEDC2
<b>Immunogen:</b>	Fusion protein of human CUEDC2
<b>Full name:</b>	CUE domain containing 2
<b>Synonyms:</b>	C10orf66; bA18I14.5
<b>SwissProt:</b>	Q9H467

<b>IHC positive control:</b>	Human gastric cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200

